



Lipika
Technologies

AUM PULSER V2.0

A PEMF Device For sound sleep and pain management



Small & Portable
Easy to Use
Zero Side effects





INTENDED USE

AUM PULSER IS USED TO ENHANCE SLEEP, PERFORMANCE, AND VITALITY

Program Modes in AUM PULSER

AUM PULSER has been programmed for the following modes for different end uses.

- **H-SLEEP:** For dealing with sleep-related problems.
- **F-SLEEP:** To fall asleep faster.
- **D-SLEEP:** For extremely sleep deprived or highly stressed.
- **AWAKE:** For working late/driving without fatigue.
- **WELLNESS:** Reduces stress level and helps to feel calmer.
- **REVITALIZE:** For pain and injury management.
- **SCHUMANN:** Connecting with Earth's natural frequency for body rejuvenation and balance.
- **MANUAL:** Set for any pulsed magnetic field frequency.

BENEFITS

AUM PULSER is designed and developed based on PEMF technology.

- PEMF is catching on as a non-invasive way to approach injuries, chronic pain, and even chronic conditions like depression and diabetes.
- This system is designed for better sleep, radically accelerated recovery, enhanced daytime performance, and incredible anti-aging effects.
- Regular use of PEMF therapy, besides promoting various healing mechanisms, has been found to have substantially beneficial neuroendocrine, neurological, and psychological effects. as well as having the ability to promote bone, tissue, and nerve regeneration.

General instruction to set up aum pulser

1. Plug the power adapter cable into the large jack at the back of the program controller.
2. Plug the power adaptor into the electricity outlet on the wall(100-240 volts AC). Use an International pin adaptor if required.
3. Plug the 2 coil sets into the remaining 2 jacks at the back of the program controller.
4. Press the "On/Off" button on the program controller, the display shows "AUM PULSER" and instantly boots to "H-SLEEP" mode(default).

NOTE: Your system is now ready to use and you need to select and set up the mode, Timer, and Amplitude.

5. Press "Enter" then the Mode changes from light to dark shaded and uses the "Up"/"Down" Arrow key to arrive at desired mode.
6. The once-desired mode is selected. Press the "Enter" key. The display now shows the mode with the timer set to 8:00 hours and Amplitude set to 10%.

NOTE: You now need to set time duration and Amplitude.

7. Press "Enter" then the Timer selection changes from light to dark shade and use the "Up" or "Down" Arrow to set the time duration for which you wish to use the device. Time can be set from 15 minutes to 12 hours with changes in multiples of 15 minutes.
8. Press "Enter" then the Amplitude selection changes from light to dark shade and use the "Up" or "Down" key to select the Amplitude. Amplitude can be set from 10% to 100% with changes of 10%.



9. Press "Enter" and the device begins to function. This is established by the blinking of the * in the display and the blue LED on the program controller.

NOTE: If you make a mistake anywhere in the steps define above, press the "Enter" button, and then restart from step 5.

The Electromagnetic coil is placed under the mattress/ pillow(always) based on a practical understanding of the point of pain and the control unit is used to set up the program mode and timer to initialize.

One User: The coil is placed under the neck and thigh region as shown below.

Two User: Each coil is placed under **each individual.**



COIL-SET PLACEMENT

Given beneath are coil-set placement location,

Body parts to be treated	Postur	Coil-set placement
Foot	Sitting	Below foot
Ankle	Sitting/Lying	Below foot/ankle
Lower Calf	Lying	Below lower calf
Knee	Sitting/Lying	Below/on top of the knee based on point of pain
Thigh	Sitting/Lying	Below/on top of thigh based on the point of pain
Lower/Middle/Upper Back	Lying	Below Lower/Middle/Upper back
Neck	Lying	Below neck
Shoulder	Lying	Below shoulder

PRECAUTIONS

To ensure the safety of the user, and the by stander, please pay attention to the following safety instructions.

- **PREGNANCY:** Pregnant women Should Not Use the device at all.
- **PACEMAKERS:** people with Heart Pacemakers, Deep Brain Stimulators, and other Electronic implants Should Not use the device at all.
- Place only the electromagnetic under the bedding. **DO NOT** put the control box unit under pillows or mattresses as overheating may occur.
- **NORTH POLE:** North side of the electromagnetic is positioned towards the body at all times.
- **SOUTH POLE:** The south pole is located where the cord ties into the bottom side of the electromagnet.
- **CLEANING:** AUM PULSER device can be gently cleaned with a piece of cloth.
- Use adequate padding, i.e. a folded towel or small pillow between body parts and coil. The coil should NOT directly touch the body parts.

DISCLAIMER

- Lipika Technologies makes no medical claims, real or implied, as to the benefit of our device and methods.
- AUM PULSER is strictly a performance Enhancement Tool utilizing Natural Field supplementation to achieve deeper states of sleep, improved performance, and recovery from pain.
- Our product is not intended to be used to diagnose, treat, cure, or prevent any disease.
- Any information provided by us is not to be construed as medical advice. The information herein has not been evaluated by FDA.
- Worldwide, there are no governmental health agencies that recognize a need to supplement natural coil fields.